



Handheld Ultrasound
Scanner Manufacturer



WIRELESS PROBE ULTRASOUND

AFFORDABLE
VISIBLE - FAST - EASY



GUANGZHOU SONOHEALTH MEDICAL
TECHNOLOGIES CO., LTD.

Address: Room 601, Building 3, No.7 Ruitai Road, Huangpu District, Guangzhou, China

E-mail: hello@sono-health.com

Telephone: +86-02-87831007

Phone/WhatsApp: +86-15013156300

Website: www.sono-health.com

YOU CAN USE IT ANYTIME & ANYWHERE

APPLICATION

Application Sites: Vascular access, Nerve

Application Directions:

- Anesthesia, ICU, intervention, tumor, catheter, emergency: Visual guidance of cvc, picc, port, long and midline catheters, radial artery, femoral artery and vein, nerve block and other puncture visual guidance.
- Renal and Hemodialysis center: Establishment and maintenance of vascular access, visual guidance of Dialysis puncture.
- Pain, rehabilitation, Physiotherapy: Visual guidance of diagnosis and treatment, such as soft tissue lysis, effusion extraction, etc.
- Plastic surgery: Visual guidance of cosmetic surgery.



FEATURES

- Integrated palm color ultrasound, built-in probe, cable aseptic problem solved.
- Integrated palm color ultrasound, built-in touch and display screen, can working independently not need connecting to tablet.
- The technique of puncture enhancement makes the puncture more reliable.
- Full sealing design: It can be short soaked and disinfected to strengthen infection control.

SPECIFICATION

- Scan Angle/Length: 30mm
- Frequency: 7/10 MHz
- Scan Mode: B, Color, PW, Volume
- Display Depth: 20~60mm
- Power: Built-in lithium battery, USB charge 3h, work 3h
- Working system: iOS, Android
- Wifi Type: 802.11n/2.4G/5G
- Dimension: 150x53x22mm
- Weight: 185g



WATERPROOF

Reduce medical accidents caused by cross infection, protect doctors and patients, and reduce hospital related expenses.



WIRELESS CHARGE

Easy and connect with a wireless power bank.



WiFi CONNECT

WiFi connect with your Android phone, iPhone or iPad.



EASY TO USE

Simple operation, short learning curve for people without ultrasound experience.



HIGH QUALITY IMAGES

See The Difference Great Image Quality Makes

